

ABX Minocal

LOT

CX 431

CAL



(Exp.) 2019-03-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.5			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.46			4.50					± 0.06
HB HGB	g/dl	13.5	13.1			13.5					± 0.2
	g/l	135	131			135					± 2
	mmol/l	8.38	8.15			8.38					± 0.12
HT HCT	%	38.1	37.7			37.8					± 1.0
	l/l	0.381	0.377			0.378					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	243			246					± 10
VPM MPV	µm ³ ; fl	8.8	8.5			8.8					± 0.5

ABX Minocal

LOT

CX 431
Rev 1

CAL

 (Exp.) **2019-03-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.5	8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.51							± 0.06
HB HGB	g/dl	13.2	13.1							± 0.2
	g/l	132	131							± 2
	mmol/l	8.20	8.14							± 0.12
HT HCT	%	37.7	37.8							± 1.0
	l/l	0.377	0.378							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	224	246							± 10
VPM MPV	µm ³ ; fl	N/A	8.9							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881