

# ABX Minocal

**LOT**

**CX 431**

**CAL**



**(Exp.) 2019-03-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.9</b>	<b>8.5</b>			<b>8.9</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.54</b>	<b>4.46</b>			<b>4.50</b>					± 0.06
HB HGB	g/dl	<b>13.5</b>	<b>13.1</b>			<b>13.5</b>					± 0.2
	g/l	<b>135</b>	<b>131</b>			<b>135</b>					± 2
	mmol/l	<b>8.38</b>	<b>8.15</b>			<b>8.38</b>					± 0.12
HT HCT	%	<b>38.1</b>	<b>37.7</b>			<b>37.8</b>					± 1.0
	l/l	<b>0.381</b>	<b>0.377</b>			<b>0.378</b>					± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>240</b>	<b>243</b>			<b>246</b>					± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>8.8</b>	<b>8.5</b>			<b>8.8</b>					± 0.5

# ABX Minocal

**LOT**

**CX 431**  
Rev 1

**CAL**

 (Exp.) **2019-03-05**  
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.5</b>	<b>8.6</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.49</b>	<b>4.51</b>							± 0.06
HB HGB	g/dl	<b>13.2</b>	<b>13.1</b>							± 0.2
	g/l	<b>132</b>	<b>131</b>							± 2
	mmol/l	<b>8.20</b>	<b>8.14</b>							± 0.12
HT HCT	%	<b>37.7</b>	<b>37.8</b>							± 1.0
	l/l	<b>0.377</b>	<b>0.378</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>224</b>	<b>246</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>8.9</b>							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881