

ABX Minocal

LOT

CX 421

CAL



(Exp.) 2018-05-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.6	9.0			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.82	4.78			4.80					± 0.06
HB HGB	g/dl	14.0	13.6			14.0					± 0.2
	g/l	140	136			140					± 2
	mmol/l	8.69	8.46			8.69					± 0.12
HT HCT	%	39.5	39.4			39.6					± 1.0
	l/l	0.395	0.394			0.396					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	243			240					± 10
VPM MPV	µm ³ ; fl	9.0	8.6			8.9					± 0.5

ABX Minocal

LOT

CX 421
Rev 1

CAL

 (Exp.) **2018-05-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.5	8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.82	4.83							± 0.06
HB HGB	g/dl	13.6	13.6							± 0.2
	g/l	136	136							± 2
	mmol/l	8.45	8.45							± 0.12
HT HCT	%	39.7	39.6							± 1.0
	l/l	0.397	0.396							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	237							± 10
VPM MPV	µm ³ ; fl	N/A	9.0							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A