

ABX Minocal

LOT

CX 450

CAL



(Exp.) 2020-10-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.8			8.7					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.46			4.44					± 0.06
HB HGB	g/dl	14.0	13.5			13.8					± 0.2
	g/l	140	135			138					± 2
	mmol/l	8.69	8.40			8.57					± 0.12
HT HCT	%	39.6	39.2			39.5					± 1.0
	l/l	0.396	0.392			0.395					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	230			235					± 10
VPM MPV	µm ³ ; fl	8.7	8.5			8.6					± 0.5

ABX Minocal

LOT

CX 450

CAL



(Exp.) 2020-10-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.5	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.43	4.46							± 0.06
HB HGB	g/dl	13.6	13.6							± 0.2
	g/l	136	136							± 2
	mmol/l	8.45	8.45							± 0.12
HT HCT	%	38.8	39.2							± 1.0
	l/l	0.388	0.392							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	240							± 10
VPM MPV	µm ³ ; fl	N/A	8.9							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881