

ABX Minocal

LOT		CX 430		CAL		(Exp.) 2019-02-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA		
		ES60	CRP 200			60 8P 60 18P						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.9			9.0					± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.63	4.58			4.57					± 0.06	
HB HGB	g/dl	13.7	13.4			13.7					± 0.2	
	g/l	137	134			137					± 2	
	mmol/l	8.51	8.34			8.51					± 0.12	
HT HCT	%	38.8	38.7			38.4					± 1.0	
	l/l	0.388	0.387			0.384					± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	234	248			247					± 10	
VPM MPV	µm ³ ; fl	8.8	8.4			8.6					± 0.5	

ABX Minocal


 (Exp.) **2019-02-05**
 (YYYY - MM - DD)

LOT		CX 430		CAL		(Exp.) 2019-02-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.7	8.8								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.60								± 0.06
HB HGB	g/dl	13.4	13.4								± 0.2
	g/l	134	134								± 2
	mmol/l	8.32	8.32								± 0.12
HT HCT	%	38.6	38.6								± 1.0
	l/l	0.386	0.386								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	227	244								± 10
VPM MPV	µm ³ ; fl	N/A	8.9								± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A