

ABX Minocal

LOT		CX 400		CAL		(Exp.) 2016-08-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
			CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.0			8.8					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.65			4.66					± 0.06
HB HGB	g/dl		13.5			13.3					± 0.2
	g/l		135			133					± 2
	mmol/l		8.40			8.25					± 0.12
HT HCT	%		39.1			39.1					± 1.0
	l/l		0.391			0.391					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		246			243					± 10
VPM MPV	µm ³ ; fl		8.4			7.9					± 0.5

ABX Minocal


 (Exp.) **2016-08-05**
 (YYYY - MM - DD)

LOT		CX 400		CAL		(Exp.)2016-08-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.7								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.70	4.68								± 0.06
HB HGB	g/dl	13.8	13.6								± 0.2
	g/l	138	136								± 2
	mmol/l	8.57	8.45								± 0.12
HT HCT	%	39.2	39.1								± 1.0
	l/l	0.392	0.391								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	233	247								± 10
VPM MPV	µm ³ ;fl	N/A	8.5								± 0.5

Ref: TEMP-0387 Rev.42 FRONT / RECTO 1300019741-A