

ABX Minocal

LOT

CX 409

CAL



(Exp.) 2017-05-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.1			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.59	4.56			4.55					± 0.06
HB HGB	g/dl	13.7	13.4			13.7					± 0.2
	g/l	137	134			137					± 2
	mmol/l	8.51	8.34			8.51					± 0.12
HT HCT	%	38.6	39.0			38.7					± 1.0
	l/l	0.386	0.390			0.387					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	254			250					± 10
VPM MPV	µm ³ ; fl	8.2	8.5			8.3					± 0.5

ABX Minocal

LOT

CX 409

CAL



(Exp.) 2017-05-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.59	4.62							± 0.06
HB HGB	g/dl	13.5	13.6							± 0.2
	g/l	135	136							± 2
	mmol/l	8.38	8.45							± 0.12
HT HCT	%	38.8	39.0							± 1.0
	l/l	0.388	0.390							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	238	250							± 10
VPM MPV	µm ³ ; fl	N/A	8.3							± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1300032436-A