

ABX Minocal

LOT

CX 428

CAL



(Exp.) 2018-12-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.3			9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40	4.43			4.39					± 0.06
HB HGB	g/dl	13.5	13.3			13.6					± 0.2
	g/l	135	133			136					± 2
	mmol/l	8.38	8.28			8.45					± 0.12
HT HCT	%	37.8	38.3			38.0					± 1.0
	l/l	0.378	0.383			0.380					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	244			250					± 10
VPM MPV	µm ³ ; fl	8.7	8.4			8.5					± 0.5

ABX Minocal

LOT

CX 428

Rev 1

CAL



(Exp.) **2018-12-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.41							± 0.06
HB HGB	g/dl	13.3	13.3							± 0.2
	g/l	133	133							± 2
	mmol/l	8.26	8.26							± 0.12
HT HCT	%	38.0	38.1							± 1.0
	l/l	0.380	0.381							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	234	245							± 10
VPM MPV	µm ³ ; fl	N/A	8.7							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A