

ABX Minocal

LOT
CX 418
CAL

(Exp.) 2018-02-05

 Rev **1**

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.4			9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.50			4.46					± 0.06
HB HGB	g/dl	14.0	13.6			13.9					± 0.2
	g/l	140	136			139					± 2
	mmol/l	8.69	8.46			8.63					± 0.12
HT HCT	%	39.5	39.6			39.2					± 1.0
	l/l	0.395	0.396			0.392					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	256	256			263					± 10
VPM MPV	µm ³ ; fl	8.5	8.1			8.4					± 0.5

ABX Minocal

LOT
CX 418
CAL

(Exp.) 2018-02-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	9.1								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.51								± 0.06
HB HGB	g/dl	13.6	13.7								± 0.2
	g/l	136	137								± 2
	mmol/l	8.45	8.51								± 0.12
HT HCT	%	39.4	39.5								± 1.0
	l/l	0.394	0.395								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	250								± 10
VPM MPV	µm ³ ; fl	N/A	8.4								± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1300032436-A