

ABX Minocal

LOT

CX 408

CAL



(Exp.) 2017-04-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4	9.4			9.4					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.55			4.47					± 0.06
HB HGB	g/dl	13.5	13.3			13.5					± 0.2
	g/l	135	133			135					± 2
	mmol/l	8.38	8.28			8.38					± 0.12
HT HCT	%	38.6	38.9			38.0					± 1.0
	l/l	0.386	0.389			0.380					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	255	257			258					± 10
VPM MPV	µm ³ ; fl	8.9	8.9			8.8					± 0.5

ABX Minocal

LOT

CX 408

CAL



(Exp.) 2017-04-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	9.1	9.1							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.53							± 0.06
HB HGB	g/dl	13.3	13.3							± 0.2
	g/l	133	133							± 2
	mmol/l	8.26	8.26							± 0.12
HT HCT	%	38.3	38.3							± 1.0
	l/l	0.383	0.383							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	246	258							± 10
VPM MPV	µm ³ ; fl	N/A	8.8							± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1300032436-A