

ABX Minocal

LOT

CX 417

CAL



(Exp.) 2018-01-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.3			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.65	4.63			4.65					± 0.06
HB HGB	g/dl	13.7	13.4			13.8					± 0.2
	g/l	137	134			138					± 2
	mmol/l	8.51	8.34			8.57					± 0.12
HT HCT	%	38.8	38.7			38.9					± 1.0
	l/l	0.388	0.387			0.389					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	255	251			260					± 10
VPM MPV	µm ³ ; fl	8.8	8.6			8.7					± 0.5

ABX Minocal

LOT

CX 417
Rev 1

CAL

 (Exp.) **2018-01-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	9.2							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.66	4.70							± 0.06
HB HGB	g/dl	13.4	13.5							± 0.2
	g/l	134	135							± 2
	mmol/l	8.32	8.38							± 0.12
HT HCT	%	38.9	39.0							± 1.0
	l/l	0.389	0.390							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	252							± 10
VPM MPV	µm ³ ; fl	N/A	8.6							± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1300032436-A