

ABX Minocal

LOT

CX 407

CAL



(Exp.) 2017-03-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.3			9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.58			4.56					± 0.06
HB HGB	g/dl	13.6	13.5			13.8					± 0.2
	g/l	136	135			138					± 2
	mmol/l	8.45	8.40			8.57					± 0.12
HT HCT	%	38.4	39.0			38.5					± 1.0
	l/l	0.384	0.390			0.385					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	243	245			248					± 10
VPM MPV	µm ³ ; fl	8.5	8.7			8.6					± 0.5

ABX Minocal

LOT

CX 407

Rev 1

CAL



(Exp.) **2017-03-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	9.0	9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.62							± 0.06
HB HGB	g/dl	13.6	13.5							± 0.2
	g/l	136	135							± 2
	mmol/l	8.45	8.38							± 0.12
HT HCT	%	38.9	38.9							± 1.0
	l/l	0.389	0.389							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	250							± 10
VPM MPV	µm ³ ; fl	N/A	8.6							± 0.5

Ref: TEMP-0387 Rev.43 FRONT / RECTO 1300027265-A