

ABX Minocal

LOT

CX 406

CAL



(Exp.) 2017-02-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.1			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.66	4.62			4.68					± 0.06
HB HGB	g/dl	13.7	13.4			13.8					± 0.2
	g/l	137	134			138					± 2
	mmol/l	8.51	8.34			8.57					± 0.12
HT HCT	%	39.1	38.8			39.3					± 1.0
	l/l	0.391	0.388			0.393					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	234	239			238					± 10
VPM MPV	µm ³ ; fl	8.8	9.0			9.2					± 0.5

ABX Minocal

LOT

CX 406

Rev 1

CAL



(Exp.) **2017-02-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.66	4.67							± 0.06
HB HGB	g/dl	13.6	13.5							± 0.2
	g/l	136	135							± 2
	mmol/l	8.45	8.38							± 0.12
HT HCT	%	39.1	39.1							± 1.0
	l/l	0.391	0.391							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	228	238							± 10
VPM MPV	µm ³ ; fl	N/A	8.8							± 0.5

Ref: TEMP-0387 Rev.43 FRONT / RECTO 1300027265-A