

ABX Minocal

LOT
CX 066

 Rev **1**
CAL

(Exp.) 2016-07-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
			CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.1			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.52			4.53					± 0.06
HB HGB	g/dl		13.4			13.4					± 0.2
	g/l		134			134					± 2
	mmol/l		8.34			8.31					± 0.12
HT HCT	%		38.4			38.6					± 1.0
	l/l		0.384			0.386					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		245			250					± 10
VPM MPV	µm ³ ; fl		8.2			8.2					± 0.5

ABX Minocal



LOT		CX 066		CAL		(Exp.) 2016-07-05 (YYYY - MM - DD)						
		Rev 1	LyseVet									
PARAMETRES PARAMETERS		UNITES UNITS		MICROS	MICROS							TOLERANCES TOLERANCE
				Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC		10 ³ /mm ³ ;10 ⁹ /l	8.6	8.8								± 0.2
GR RBC		10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.58								± 0.06
HB HGB		g/dl	13.6	13.5								± 0.2
		g/l	136	135								± 2
		mmol/l	8.45	8.38								± 0.12
HT HCT		%	38.7	38.7								± 1.0
		l/l	0.387	0.387								± 0.010
PLAQ. PLTS		10 ³ /mm ³ ;10 ⁹ /l	237	253								± 10
VPM MPV		µm ³ ; fl	N/A	8.4								± 0.5

Ref: TEMP-0387 Rev.42 FRONT / RECTO 1300019741-A