

ABX Minocal

LOT

CX 471

CAL



(Exp.) 2022-07-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.8			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.49			4.49					± 0.06
HB HGB	g/dl	13.5	13.3			13.5					± 0.2
	g/l	135	133			135					± 2
	mmol/l	8.38	8.28			8.38					± 0.12
HT HCT	%	38.3	38.6			39.1					± 1.0
	l/l	0.383	0.386			0.391					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	244			250					± 10
VPM MPV	µm ³ ; fl	8.8	8.6			9.0					± 0.5

ABX Minocal

LOT

CX 471

CAL



(Exp.) 2022-07-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.6	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.48							± 0.06
HB HGB	g/dl	13.2	13.2							± 0.2
	g/l	132	132							± 2
	mmol/l	8.20	8.20							± 0.12
HT HCT	%	37.7	38.1							± 1.0
	l/l	0.377	0.381							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	238	249							± 10
VPM MPV	µm ³ ; fl	N/A	9.3							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881