

# ABX Minocal



**LOT** CX 061  
Rev 1

**CAL**

(Exp.) 2011-07-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE	
		MICROS	MICROS										
		Scil Vet ABC	Scil Vet ABC Plus ESV 60										
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	<b>8.8</b>	<b>8.8</b>										± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l	<b>4.62</b>	<b>4.65</b>										± 0.06
HB HGB	g/dl	<b>13.6</b>	<b>13.4</b>										± 0.2
	g/l	<b>136</b>	<b>134</b>										± 2
	mmol/l	<b>8.45</b>	<b>8.34</b>										± 0.12
HT HCT	%	<b>39.9</b>	<b>39.6</b>										± 1.0
	l/l	<b>0.399</b>	<b>0.396</b>										± 0.010
PLAQ. PLTS	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	<b>244</b>	<b>240</b>										± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>8.2</b>										± 0.5

Ref: FORM-0387 Rev.17 FRONT / RECTO 9930044-B