

ABX Minocal

LOT

CX 463

CAL



(Exp.) 2021-11-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.8			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.37	4.36			4.39					± 0.06
HB HGB	g/dl	13.3	13.0			13.4					± 0.2
	g/l	133	130			134					± 2
	mmol/l	8.26	8.09			8.32					± 0.12
HT HCT	%	37.4	37.5			38.4					± 1.0
	l/l	0.374	0.375			0.384					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	232	245			245					± 10
VPM MPV	µm ³ ; fl	7.8	7.6			7.9					± 0.5

ABX Minocal

LOT

CX 463

CAL



(Exp.) 2021-11-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.6	8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.34	4.33							± 0.06
HB HGB	g/dl	12.9	12.9							± 0.2
	g/l	129	129							± 2
	mmol/l	8.01	8.01							± 0.12
HT HCT	%	36.9	36.8							± 1.0
	l/l	0.369	0.368							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	241							± 10
VPM MPV	µm ³ ; fl	N/A	8.2							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881