

# ABX Minocal

**LOT**

**CX 463**

**CAL**



**(Exp.) 2021-11-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.8</b>	<b>8.8</b>			<b>9.0</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.37</b>	<b>4.36</b>			<b>4.39</b>					± 0.06
HB HGB	g/dl	<b>13.3</b>	<b>13.0</b>			<b>13.4</b>					± 0.2
	g/l	<b>133</b>	<b>130</b>			<b>134</b>					± 2
	mmol/l	<b>8.26</b>	<b>8.09</b>			<b>8.32</b>					± 0.12
HT HCT	%	<b>37.4</b>	<b>37.5</b>			<b>38.4</b>					± 1.0
	l/l	<b>0.374</b>	<b>0.375</b>			<b>0.384</b>					± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>232</b>	<b>245</b>			<b>245</b>					± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>7.8</b>	<b>7.6</b>			<b>7.9</b>					± 0.5

# ABX Minocal

**LOT**

**CX 463**

**CAL**



**(Exp.) 2021-11-05**

(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.6</b>	<b>8.6</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.34</b>	<b>4.33</b>							± 0.06
HB HGB	g/dl	<b>12.9</b>	<b>12.9</b>							± 0.2
	g/l	<b>129</b>	<b>129</b>							± 2
	mmol/l	<b>8.01</b>	<b>8.01</b>							± 0.12
HT HCT	%	<b>36.9</b>	<b>36.8</b>							± 1.0
	l/l	<b>0.369</b>	<b>0.368</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>230</b>	<b>241</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>8.2</b>							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881