

ABX Minocal

LOT

CX 453

CAL



(Exp.) 2021-01-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.0			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.55			4.50					± 0.06
HB HGB	g/dl	14.1	13.9			14.1					± 0.2
	g/l	141	139			141					± 2
	mmol/l	8.76	8.64			8.76					± 0.12
HT HCT	%	39.9	40.0			40.1					± 1.0
	l/l	0.399	0.400			0.401					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	232	240			240					± 10
VPM MPV	µm ³ ; fl	8.7	8.2			8.9					± 0.5

ABX Minocal

LOT

CX 453

CAL



(Exp.) 2021-01-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.50							± 0.06
HB HGB	g/dl	13.9	13.8							± 0.2
	g/l	139	138							± 2
	mmol/l	8.63	8.57							± 0.12
HT HCT	%	40.2	40.0							± 1.0
	l/l	0.402	0.400							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	238	244							± 10
VPM MPV	µm ³ ; fl	N/A	9.0							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881