

ABX Minocal

LOT

CX 443

CAL



(Exp.) 2020-03-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	8.9			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.50			4.54					± 0.06
HB HGB	g/dl	14.1	13.6			14.1					± 0.2
	g/l	141	136			141					± 2
	mmol/l	8.76	8.46			8.76					± 0.12
HT HCT	%	39.6	39.4			39.7					± 1.0
	l/l	0.396	0.394			0.397					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	250	235			245					± 10
VPM MPV	µm ³ ; fl	8.3	8.1			8.8					± 0.5

ABX Minocal

LOT

CX 443

CAL



(Exp.) 2020-03-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.51							± 0.06
HB HGB	g/dl	13.8	13.7							± 0.2
	g/l	138	137							± 2
	mmol/l	8.57	8.51							± 0.12
HT HCT	%	39.5	39.2							± 1.0
	l/l	0.395	0.392							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	233	252							± 10
VPM MPV	µm ³ ; fl	N/A	8.8							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881