

ABX Minocal

LOT

CX 433

CAL



(Exp.) 2019-05-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.8			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.38			4.40					± 0.06
HB HGB	g/dl	13.4	13.2			13.4					± 0.2
	g/l	134	132			134					± 2
	mmol/l	8.32	8.22			8.32					± 0.12
HT HCT	%	37.8	37.7			37.7					± 1.0
	l/l	0.378	0.377			0.377					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	234	248			249					± 10
VPM MPV	µm ³ ; fl	8.7	8.0			8.5					± 0.5

ABX Minocal

LOT

CX 433

Rev 1

CAL



(Exp.) **2019-05-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.6	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40	4.38							± 0.06
HB HGB	g/dl	13.2	13.1							± 0.2
	g/l	132	131							± 2
	mmol/l	8.20	8.14							± 0.12
HT HCT	%	37.8	37.7							± 1.0
	l/l	0.378	0.377							± 0.010
PLAQ. PLTS	10 ⁹ /mm ³ ;10 ⁹ /l	232	242							± 10
VPM MPV	µm ³ ; fl	N/A	8.6							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881