


ABX Minocal

LOT CX 413
Rev 2

CAL

 (Exp.) 2017-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.1			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.69	4.68			4.64					± 0.06
HB HGB	g/dl	14.0	13.5			13.9					± 0.2
	g/l	140	135			139					± 2
	mmol/l	8.69	8.40			8.63					± 0.12
HT HCT	%	39.6	39.8			39.4					± 1.0
	l/l	0.396	0.398			0.394					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	254	252			256					± 10
VPM MPV	µm ³ ; fl	8.9	8.4			8.7					± 0.5

ABX Minocal



LOT CX 413
Rev 2

CAL

(Exp.) 2017-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	8.8	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.67	4.67							± 0.06
HB HGB	g/dl	13.6	13.7							± 0.2
	g/l	136	137							± 2
	mmol/l	8.45	8.51							± 0.12
HT HCT	%	39.2	39.2							± 1.0
	l/l	0.392	0.392							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	241	252							± 10
VPM MPV	µm ³ ; fl	N/A	8.6							± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1300032436-A