

ABX Minocal

LOT

CX 462

CAL



(Exp.) 2021-10-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.0			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.60			4.56					± 0.06
HB HGB	g/dl	13.7	13.5			13.7					± 0.2
	g/l	137	135			137					± 2
	mmol/l	8.51	8.40			8.51					± 0.12
HT HCT	%	39.0	39.3			39.7					± 1.0
	l/l	0.390	0.393			0.397					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	236	250			245					± 10
VPM MPV	µm ³ ; fl	9.0	8.7			9.2					± 0.5

ABX Minocal

LOT

CX 462

CAL



(Exp.) 2021-10-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.52							± 0.06
HB HGB	g/dl	13.3	13.4							± 0.2
	g/l	133	134							± 2
	mmol/l	8.26	8.32							± 0.12
HT HCT	%	38.6	38.6							± 1.0
	l/l	0.386	0.386							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	244							± 10
VPM MPV	µm ³ ; fl	N/A	9.3							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881