

ABX Minocal

LOT		CX 452		CAL		(Exp.) 2020-12-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA		
		ES60	CRP 200			60 8P 60 18P						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3	9.2			9.2					± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.55			4.55					± 0.06	
HB HGB	g/dl	14.0	13.8			13.9					± 0.2	
	g/l	140	138			139					± 2	
	mmol/l	8.69	8.58			8.63					± 0.12	
HT HCT	%	39.4	40.0			40.0					± 1.0	
	l/l	0.394	0.400			0.400					± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	237			240					± 10	
VPM MPV	µm ³ ; fl	8.3	7.9			8.6					± 0.5	

ABX Minocal



LOT		CX 452		CAL		(Exp.)2020-12-05 (YYYY - MM - DD)									
		Rev 1	LyseVet												
PARAMETRES PARAMETERS		UNITES UNITS										TOLERANCES TOLERANCE			
				MICROS		MICROS									
				Scil Vet ABC	Scil Vet ABC Plus ESV 60										
GB WBC		10 ³ /mm ³ ;10 ⁹ /l	8.6	8.8								± 0.2			
GR RBC		10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.50								± 0.06			
HB HGB		g/dl	13.7	13.7								± 0.2			
		g/l	137	137								± 2			
		mmol/l	8.51	8.51								± 0.12			
HT HCT		%	39.5	39.6								± 1.0			
		l/l	0.395	0.396								± 0.010			
PLAQ. PLTS		10 ³ /mm ³ ;10 ⁹ /l	230	238								± 10			
VPM MPV		µm ³ ; fl	N/A	8.7								± 0.5			

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881