

ABX Minocal

LOT

CX 432

CAL



(Exp.) 2019-04-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	8.6			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.38	4.28			4.34					± 0.06
HB HGB	g/dl	13.4	13.0			13.3					± 0.2
	g/l	134	130			133					± 2
	mmol/l	8.32	8.09			8.26					± 0.12
HT HCT	%	37.7	37.2			37.3					± 1.0
	l/l	0.377	0.372			0.373					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	236	245			249					± 10
VPM MPV	µm ³ ; fl	9.0	8.6			8.8					± 0.5

ABX Minocal

LOT

CX 432

Rev 1

CAL



(Exp.) **2019-04-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.35	4.36							± 0.06
HB HGB	g/dl	13.2	13.0							± 0.2
	g/l	132	130							± 2
	mmol/l	8.20	8.07							± 0.12
HT HCT	%	37.6	37.5							± 1.0
	l/l	0.376	0.375							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	250							± 10
VPM MPV	µm ³ ; fl	N/A	9.1							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881