

ABX Minocal

LOT

CX 412

CAL



(Exp.) 2017-08-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3	9.3			9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.55			4.50					± 0.06
HB HGB	g/dl	14.0	13.5			13.9					± 0.2
	g/l	140	135			139					± 2
	mmol/l	8.69	8.40			8.63					± 0.12
HT HCT	%	39.4	39.4			39.2					± 1.0
	l/l	0.394	0.394			0.392					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	232			237					± 10
VPM MPV	µm ³ ; fl	8.4	8.1			8.1					± 0.5

ABX Minocal

LOT

CX 412

Rev 1

CAL



(Exp.) 2017-08-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.54							± 0.06
HB HGB	g/dl	13.6	13.6							± 0.2
	g/l	136	136							± 2
	mmol/l	8.45	8.45							± 0.12
HT HCT	%	39.2	39.0							± 1.0
	l/l	0.392	0.390							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	225	234							± 10
VPM MPV	µm ³ ; fl	N/A	8.0							± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1300032436-A