

ABX Minocal

LOT

CX 461

CAL



(Exp.) 2021-09-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.6			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.60			4.58					± 0.06
HB HGB	g/dl	14.2	13.8			14.3					± 0.2
	g/l	142	138			143					± 2
	mmol/l	8.82	8.58			8.88					± 0.12
HT HCT	%	39.8	40.0			41.0					± 1.0
	l/l	0.398	0.400			0.410					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	228	235			236					± 10
VPM MPV	µm ³ ; fl	8.0	7.7			8.0					± 0.5

ABX Minocal

LOT

CX 461

CAL



(Exp.) 2021-09-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.6	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.55							± 0.06
HB HGB	g/dl	13.7	13.9							± 0.2
	g/l	137	139							± 2
	mmol/l	8.51	8.63							± 0.12
HT HCT	%	39.1	39.6							± 1.0
	l/l	0.391	0.396							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	225	231							± 10
VPM MPV	µm ³ ; fl	N/A	8.3							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881