

ABX Minocal

LOT

CX 451

CAL



(Exp.) 2020-11-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3	9.2			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.47			4.50					± 0.06
HB HGB	g/dl	14.0	13.7			13.9					± 0.2
	g/l	140	137			139					± 2
	mmol/l	8.69	8.52			8.63					± 0.12
HT HCT	%	39.6	39.6			39.6					± 1.0
	l/l	0.396	0.396			0.396					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	236	242			245					± 10
VPM MPV	µm ³ ; fl	8.4	8.0			8.5					± 0.5

ABX Minocal

LOT

CX 451

CAL



(Exp.) 2020-11-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.48							± 0.06
HB HGB	g/dl	13.7	13.7							± 0.2
	g/l	137	137							± 2
	mmol/l	8.51	8.51							± 0.12
HT HCT	%	39.4	39.6							± 1.0
	l/l	0.394	0.396							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	234	247							± 10
VPM MPV	µm ³ ; fl	N/A	8.7							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881