

ABX Minocal

LOT		CX 421		CAL		(Exp.) 2018-05-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.6	9.0			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.82	4.78			4.80					± 0.06
HB HGB	g/dl	14.0	13.6			14.0					± 0.2
	g/l	140	136			140					± 2
	mmol/l	8.69	8.46			8.69					± 0.12
HT HCT	%	39.5	39.4			39.6					± 1.0
	l/l	0.395	0.394			0.396					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	243			240					± 10
VPM MPV	µm ³ ; fl	9.0	8.6			8.9					± 0.5

ABX Minocal

LOT		CX 421		CAL		(Exp.)2018-05-05 (YYYY - MM - DD)							
		Rev 1		LyseVet									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS								TOLERANCES TOLERANCE	
		Scil Vet ABC	Scil Vet ABC Plus ESV 60										
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.5	8.6								± 0.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.82	4.83								± 0.06		
HB HGB	g/dl	13.6	13.6								± 0.2		
	g/l	136	136								± 2		
	mmol/l	8.45	8.45								± 0.12		
HT HCT	%	39.7	39.6								± 1.0		
	l/l	0.397	0.396								± 0.010		
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	237								± 10		
VPM MPV	µm ³ ; fl	N/A	9.0								± 0.5		

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A