

# ABX Minocal

**LOT**

**CX 421**

**CAL**



**(Exp.) 2018-05-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.6</b>	<b>9.0</b>			<b>8.9</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.82</b>	<b>4.78</b>			<b>4.80</b>					± 0.06
HB HGB	g/dl	<b>14.0</b>	<b>13.6</b>			<b>14.0</b>					± 0.2
	g/l	<b>140</b>	<b>136</b>			<b>140</b>					± 2
	mmol/l	<b>8.69</b>	<b>8.46</b>			<b>8.69</b>					± 0.12
HT HCT	%	<b>39.5</b>	<b>39.4</b>			<b>39.6</b>					± 1.0
	l/l	<b>0.395</b>	<b>0.394</b>			<b>0.396</b>					± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>245</b>	<b>243</b>			<b>240</b>					± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>9.0</b>	<b>8.6</b>			<b>8.9</b>					± 0.5

# ABX Minocal

**LOT**

**CX 421**  
Rev 1

**CAL**

 (Exp.) **2018-05-05**  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.5</b>	<b>8.6</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.82</b>	<b>4.83</b>							± 0.06
HB HGB	g/dl	<b>13.6</b>	<b>13.6</b>							± 0.2
	g/l	<b>136</b>	<b>136</b>							± 2
	mmol/l	<b>8.45</b>	<b>8.45</b>							± 0.12
HT HCT	%	<b>39.7</b>	<b>39.6</b>							± 1.0
	l/l	<b>0.397</b>	<b>0.396</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>230</b>	<b>237</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>9.0</b>							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A