

ABX Minocal

LOT

CX 411

CAL



(Exp.) 2017-07-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.0			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.59	4.57			4.59					± 0.06
HB HGB	g/dl	13.9	13.4			13.9					± 0.2
	g/l	139	134			139					± 2
	mmol/l	8.63	8.34			8.63					± 0.12
HT HCT	%	39.5	39.1			39.0					± 1.0
	l/l	0.395	0.391			0.390					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	247	250			253					± 10
VPM MPV	µm ³ ; fl	8.5	8.2			8.2					± 0.5

ABX Minocal

LOT

CX 411
Rev 1

CAL



(Exp.) 2017-07-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.61	4.61							± 0.06
HB HGB	g/dl	13.6	13.6							± 0.2
	g/l	136	136							± 2
	mmol/l	8.45	8.45							± 0.12
HT HCT	%	39.2	39.0							± 1.0
	l/l	0.392	0.390							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	250							± 10
VPM MPV	µm ³ ; fl	N/A	8.4							± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1300032436-A