

ABX Minocal

LOT

CX 440

CAL



(Exp.) 2019-12-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	8.5			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.45			4.49					± 0.06
HB HGB	g/dl	13.8	13.5			13.8					± 0.2
	g/l	138	135			138					± 2
	mmol/l	8.57	8.40			8.57					± 0.12
HT HCT	%	39.1	38.7			39.0					± 1.0
	l/l	0.391	0.387			0.390					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	243	240			240					± 10
VPM MPV	µm ³ ; fl	9.0	8.6			9.3					± 0.5

ABX Minocal

LOT

CX 440

CAL



(Exp.) 2019-12-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.50							± 0.06
HB HGB	g/dl	13.6	13.5							± 0.2
	g/l	136	135							± 2
	mmol/l	8.45	8.38							± 0.12
HT HCT	%	38.9	38.8							± 1.0
	l/l	0.389	0.388							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	250							± 10
VPM MPV	µm ³ ; fl	N/A	9.4							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881