

ABX Minocal

LOT

CX 430

CAL



(Exp.) 2019-02-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.9			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.63	4.58			4.57					± 0.06
HB HGB	g/dl	13.7	13.4			13.7					± 0.2
	g/l	137	134			137					± 2
	mmol/l	8.51	8.34			8.51					± 0.12
HT HCT	%	38.8	38.7			38.4					± 1.0
	l/l	0.388	0.387			0.384					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	234	248			247					± 10
VPM MPV	µm ³ ; fl	8.8	8.4			8.6					± 0.5

ABX Minocal

LOT

CX 430

CAL



(Exp.) 2019-02-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.60							± 0.06
HB HGB	g/dl	13.4	13.4							± 0.2
	g/l	134	134							± 2
	mmol/l	8.32	8.32							± 0.12
HT HCT	%	38.6	38.6							± 1.0
	l/l	0.386	0.386							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	227	244							± 10
VPM MPV	µm ³ ; fl	N/A	8.9							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A