

# ABX Minocal

**LOT**

**CX 429**

**CAL**



**(Exp.) 2019-01-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.0</b>	<b>8.9</b>			<b>9.1</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.60</b>	<b>4.52</b>			<b>4.52</b>					± 0.06
HB HGB	g/dl	<b>14.1</b>	<b>13.7</b>			<b>14.1</b>					± 0.2
	g/l	<b>141</b>	<b>137</b>			<b>141</b>					± 2
	mmol/l	<b>8.76</b>	<b>8.52</b>			<b>8.76</b>					± 0.12
HT HCT	%	<b>39.6</b>	<b>39.3</b>			<b>39.1</b>					± 1.0
	l/l	<b>0.396</b>	<b>0.393</b>			<b>0.391</b>					± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>225</b>	<b>232</b>			<b>233</b>					± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>8.0</b>	<b>7.6</b>			<b>8.0</b>					± 0.5

# ABX Minocal

**LOT**

**CX 429**

**CAL**



**(Exp.) 2019-01-05**  
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.9</b>	<b>8.9</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.58</b>	<b>4.60</b>							± 0.06
HB HGB	g/dl	<b>14.0</b>	<b>13.9</b>							± 0.2
	g/l	<b>140</b>	<b>139</b>							± 2
	mmol/l	<b>8.69</b>	<b>8.63</b>							± 0.12
HT HCT	%	<b>39.8</b>	<b>40.0</b>							± 1.0
	l/l	<b>0.398</b>	<b>0.400</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>220</b>	<b>238</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>8.1</b>							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A