

ABX Minocal

LOT

CX 429

CAL



(Exp.) 2019-01-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.9			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.60	4.52			4.52					± 0.06
HB HGB	g/dl	14.1	13.7			14.1					± 0.2
	g/l	141	137			141					± 2
	mmol/l	8.76	8.52			8.76					± 0.12
HT HCT	%	39.6	39.3			39.1					± 1.0
	l/l	0.396	0.393			0.391					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	225	232			233					± 10
VPM MPV	µm ³ ; fl	8.0	7.6			8.0					± 0.5

ABX Minocal

LOT

CX 429

CAL



(Exp.) 2019-01-05

(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.60							± 0.06
HB HGB	g/dl	14.0	13.9							± 0.2
	g/l	140	139							± 2
	mmol/l	8.69	8.63							± 0.12
HT HCT	%	39.8	40.0							± 1.0
	l/l	0.398	0.400							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	220	238							± 10
VPM MPV	µm ³ ; fl	N/A	8.1							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A