


ABX Minocal

LOT CX 404
Rev 1

CAL

 (Exp.) 2016-12-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.5			9.6					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.60			4.61					± 0.06
HB HGB	g/dl		13.5			13.8					± 0.2
	g/l		135			138					± 2
	mmol/l		8.40			8.57					± 0.12
HT HCT	%		39.2			39.4					± 1.0
	l/l		0.392			0.394					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		242			248					± 10
VPM MPV	µm ³ ; fl		8.4			8.4					± 0.5

ABX Minocal

LOT

CX 404

Rev 1

CAL



(Exp.) **2016-12-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	9.3	9.3							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.65	4.65							± 0.06
HB HGB	g/dl	13.7	13.6							± 0.2
	g/l	137	136							± 2
	mmol/l	8.51	8.45							± 0.12
HT HCT	%	39.5	39.5							± 1.0
	l/l	0.395	0.395							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	233	244							± 10
VPM MPV	µm ³ ; fl	N/A	8.3							± 0.5

Ref: TEMP-0387 Rev.43 FRONT / RECTO 1300027265-A