

ABX Minocal

LOT

CX 424

CAL



(Exp.) 2018-08-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.0			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.55			4.60					± 0.06
HB HGB	g/dl	13.9	13.4			13.9					± 0.2
	g/l	139	134			139					± 2
	mmol/l	8.63	8.34			8.63					± 0.12
HT HCT	%	38.7	38.7			39.0					± 1.0
	l/l	0.387	0.387			0.390					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	248			253					± 10
VPM MPV	µm ³ ; fl	8.9	8.5			8.6					± 0.5

ABX Minocal

LOT

CX 424

Rev 1

CAL



(Exp.) **2018-08-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.58							± 0.06
HB HGB	g/dl	13.5	13.5							± 0.2
	g/l	135	135							± 2
	mmol/l	8.38	8.38							± 0.12
HT HCT	%	38.8	38.8							± 1.0
	l/l	0.388	0.388							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	244							± 10
VPM MPV	µm ³ ; fl	N/A	8.8							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A