

ABX Minocal

LOT		CX 466		CAL		(Exp.) 2022-02-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA		
		ES60	CRP 200			60 8P 60 18P						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	9.0			9.0					± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.50			4.42					± 0.06	
HB HGB	g/dl	13.9	13.7			14.0					± 0.2	
	g/l	139	137			140					± 2	
	mmol/l	8.63	8.52			8.69					± 0.12	
HT HCT	%	39.2	39.6			40.2					± 1.0	
	l/l	0.392	0.396			0.402					± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	245			241					± 10	
VPM MPV	µm ³ ; fl	8.5	8.0			8.8					± 0.5	

ABX Minocal


 (Exp.) **2022-02-05**
 (YYYY - MM - DD)

LOT		CX 466		CAL		(Exp.) 2022-02-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.5	8.6								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.42								± 0.06
HB HGB	g/dl	13.5	13.5								± 0.2
	g/l	135	135								± 2
	mmol/l	8.38	8.38								± 0.12
HT HCT	%	38.7	38.5								± 1.0
	l/l	0.387	0.385								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	233	247								± 10
VPM MPV	µm ³ ; fl	N/A	9.1								± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881