

# ABX Minocal

**LOT**

**CX 470**

**CAL**



**(Exp.) 2022-06-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.1</b>	<b>9.2</b>			<b>9.3</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.42</b>	<b>4.50</b>			<b>4.46</b>					± 0.06
HB HGB	g/dl	<b>13.8</b>	<b>13.7</b>			<b>14.0</b>					± 0.2
	g/l	<b>138</b>	<b>137</b>			<b>140</b>					± 2
	mmol/l	<b>8.57</b>	<b>8.52</b>			<b>8.69</b>					± 0.12
HT HCT	%	<b>38.9</b>	<b>39.8</b>			<b>40.6</b>					± 1.0
	l/l	<b>0.389</b>	<b>0.398</b>			<b>0.406</b>					± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>232</b>	<b>241</b>			<b>240</b>					± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>9.0</b>	<b>8.4</b>			<b>9.0</b>					± 0.5

# ABX Minocal

**LOT**

**CX 470**

**CAL**



**(Exp.) 2022-06-05**  
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.8</b>	<b>9.0</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.42</b>	<b>4.46</b>							± 0.06
HB HGB	g/dl	<b>13.5</b>	<b>13.6</b>							± 0.2
	g/l	<b>135</b>	<b>136</b>							± 2
	mmol/l	<b>8.38</b>	<b>8.45</b>							± 0.12
HT HCT	%	<b>38.9</b>	<b>39.2</b>							± 1.0
	l/l	<b>0.389</b>	<b>0.392</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>237</b>	<b>245</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>9.2</b>							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881