

ABX Minocal

LOT

CX 459

CAL



(Exp.) 2021-07-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.8			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.56			4.53					± 0.06
HB HGB	g/dl	13.7	13.5			13.7					± 0.2
	g/l	137	135			137					± 2
	mmol/l	8.51	8.40			8.51					± 0.12
HT HCT	%	38.4	39.0			39.4					± 1.0
	l/l	0.384	0.390			0.394					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	254			248					± 10
VPM MPV	µm ³ ; fl	8.6	8.3			9.8					± 0.5

ABX Minocal

LOT

CX 459

CAL



(Exp.) 2021-07-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.54							± 0.06
HB HGB	g/dl	13.5	13.4							± 0.2
	g/l	135	134							± 2
	mmol/l	8.38	8.32							± 0.12
HT HCT	%	38.7	38.6							± 1.0
	l/l	0.387	0.386							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	250	255							± 10
VPM MPV	µm ³ ; fl	N/A	9.0							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881