

ABX Minocal

LOT

CX 449

CAL



(Exp.) 2020-09-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	8.8			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.50			4.42					± 0.06
HB HGB	g/dl	13.9	13.5			13.6					± 0.2
	g/l	139	135			136					± 2
	mmol/l	8.63	8.40			8.45					± 0.12
HT HCT	%	39.2	39.2			38.9					± 1.0
	l/l	0.392	0.392			0.389					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	233			233					± 10
VPM MPV	µm ³ ; fl	8.1	7.5			8.0					± 0.5

ABX Minocal

LOT

CX 449

CAL



(Exp.) 2020-09-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.48							± 0.06
HB HGB	g/dl	13.5	13.5							± 0.2
	g/l	135	135							± 2
	mmol/l	8.38	8.38							± 0.12
HT HCT	%	38.8	39.0							± 1.0
	l/l	0.388	0.390							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	234							± 10
VPM MPV	µm ³ ; fl	N/A	8.2							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881