

ABX Minocal

LOT

CX 438

CAL



(Exp.) 2019-10-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.4			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.44			4.47					± 0.06
HB HGB	g/dl	13.6	13.3			13.6					± 0.2
	g/l	136	133			136					± 2
	mmol/l	8.45	8.28			8.45					± 0.12
HT HCT	%	38.5	38.2			38.4					± 1.0
	l/l	0.385	0.382			0.384					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	240			248					± 10
VPM MPV	µm ³ ; fl	9.0	8.6			9.0					± 0.5

ABX Minocal

LOT

CX 438

CAL



(Exp.) 2019-10-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.6	8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.45							± 0.06
HB HGB	g/dl	13.5	13.4							± 0.2
	g/l	135	134							± 2
	mmol/l	8.38	8.32							± 0.12
HT HCT	%	38.7	38.5							± 1.0
	l/l	0.387	0.385							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	228	252							± 10
VPM MPV	µm ³ ; fl	N/A	9.0							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881