

ABX Minocal

LOT

CX 467

CAL



(Exp.) 2022-03-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.1			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.45	4.50			4.40					± 0.06
HB HGB	g/dl	13.5	13.2			13.4					± 0.2
	g/l	135	132			134					± 2
	mmol/l	8.38	8.22			8.32					± 0.12
HT HCT	%	37.8	38.5			38.3					± 1.0
	l/l	0.378	0.385			0.383					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	255			250					± 10
VPM MPV	µm ³ ; fl	8.2	7.8			8.8					± 0.5

ABX Minocal

LOT

CX 467

CAL



(Exp.) 2022-03-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40	4.40							± 0.06
HB HGB	g/dl	13.0	13.1							± 0.2
	g/l	130	131							± 2
	mmol/l	8.07	8.14							± 0.12
HT HCT	%	37.0	37.4							± 1.0
	l/l	0.370	0.374							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	250							± 10
VPM MPV	µm ³ ; fl	N/A	8.6							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881