

ABX Minocal

LOT

CX 457

CAL



(Exp.) 2021-05-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.9			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.47			4.56					± 0.06
HB HGB	g/dl	14.0	13.5			13.9					± 0.2
	g/l	140	135			139					± 2
	mmol/l	8.69	8.40			8.63					± 0.12
HT HCT	%	39.4	38.9			40.1					± 1.0
	l/l	0.394	0.389			0.401					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	228	239			248					± 10
VPM MPV	µm ³ ; fl	8.7	8.6			9.4					± 0.5

ABX Minocal

LOT

CX 457

CAL



(Exp.) 2021-05-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.6	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.50							± 0.06
HB HGB	g/dl	13.7	13.5							± 0.2
	g/l	137	135							± 2
	mmol/l	8.51	8.38							± 0.12
HT HCT	%	38.9	38.7							± 1.0
	l/l	0.389	0.387							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	237							± 10
VPM MPV	µm ³ ; fl	N/A	9.1							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881