

ABX Minocal

LOT

CX 437

CAL



(Exp.) 2019-09-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.5			8.8					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.48			4.47					± 0.06
HB HGB	g/dl	13.6	13.4			13.6					± 0.2
	g/l	136	134			136					± 2
	mmol/l	8.45	8.34			8.45					± 0.12
HT HCT	%	38.4	38.8			38.4					± 1.0
	l/l	0.384	0.388			0.384					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	242			245					± 10
VPM MPV	µm ³ ; fl	8.5	8.0			8.4					± 0.5

ABX Minocal

LOT

CX 437

Rev 1

CAL



(Exp.) **2019-09-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.5	8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.46							± 0.06
HB HGB	g/dl	13.5	13.3							± 0.2
	g/l	135	133							± 2
	mmol/l	8.38	8.26							± 0.12
HT HCT	%	38.7	38.4							± 1.0
	l/l	0.387	0.384							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	242							± 10
VPM MPV	µm ³ ; fl	N/A	8.5							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881