

# ABX Minocal

**LOT**

**CX 427**

**CAL**



**(Exp.) 2018-11-05**  
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.2</b>	<b>9.2</b>			<b>9.3</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.56</b>	<b>4.54</b>			<b>4.54</b>					± 0.06
HB HGB	g/dl	<b>13.8</b>	<b>13.4</b>			<b>13.7</b>					± 0.2
	g/l	<b>138</b>	<b>134</b>			<b>137</b>					± 2
	mmol/l	<b>8.57</b>	<b>8.34</b>			<b>8.51</b>					± 0.12
HT HCT	%	<b>38.8</b>	<b>38.6</b>			<b>38.6</b>					± 1.0
	l/l	<b>0.388</b>	<b>0.386</b>			<b>0.386</b>					± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>237</b>	<b>242</b>			<b>245</b>					± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>9.5</b>	<b>9.1</b>			<b>9.4</b>					± 0.5

# ABX Minocal

**LOT**

**CX 427**  
Rev 1

**CAL**



**(Exp.) 2018-11-05**  
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.0</b>	<b>9.1</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.56</b>	<b>4.56</b>							± 0.06
HB HGB	g/dl	<b>13.6</b>	<b>13.5</b>							± 0.2
	g/l	<b>136</b>	<b>135</b>							± 2
	mmol/l	<b>8.45</b>	<b>8.38</b>							± 0.12
HT HCT	%	<b>38.8</b>	<b>39.0</b>							± 1.0
	l/l	<b>0.388</b>	<b>0.390</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>228</b>	<b>235</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>9.6</b>							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A