

ABX Minocal

LOT

CX 434

CAL



(Exp.) 2019-06-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.7			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.48			4.47					± 0.06
HB HGB	g/dl	13.4	13.1			13.4					± 0.2
	g/l	134	131			134					± 2
	mmol/l	8.32	8.15			8.32					± 0.12
HT HCT	%	38.4	38.1			38.0					± 1.0
	l/l	0.384	0.381			0.380					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	240			243					± 10
VPM MPV	µm ³ ; fl	9.1	8.6			8.8					± 0.5

ABX Minocal

LOT

CX 434

CAL



(Exp.) 2019-06-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.6	8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.50							± 0.06
HB HGB	g/dl	13.2	13.1							± 0.2
	g/l	132	131							± 2
	mmol/l	8.20	8.14							± 0.12
HT HCT	%	38.0	38.3							± 1.0
	l/l	0.380	0.383							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	225	240							± 10
VPM MPV	µm ³ ; fl	N/A	9.1							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881