

ABX Minocal

LOT

CX 436

CAL



(Exp.) 2019-08-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.4			8.8					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.55			4.55					± 0.06
HB HGB	g/dl	13.8	13.5			13.7					± 0.2
	g/l	138	135			137					± 2
	mmol/l	8.57	8.40			8.51					± 0.12
HT HCT	%	38.8	38.7			38.7					± 1.0
	l/l	0.388	0.387			0.387					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	227	237			235					± 10
VPM MPV	µm ³ ; fl	8.9	8.4			9.1					± 0.5

ABX Minocal

LOT

CX 436

CAL



(Exp.) 2019-08-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.4	8.5							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.55							± 0.06
HB HGB	g/dl	13.6	13.4							± 0.2
	g/l	136	134							± 2
	mmol/l	8.45	8.32							± 0.12
HT HCT	%	38.8	38.7							± 1.0
	l/l	0.388	0.387							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	220	238							± 10
VPM MPV	µm ³ ; fl	N/A	8.9							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881