

ABX Minocal

LOT		CX 426		CAL		(Exp.) 2018-10-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.1			9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.51			4.50					± 0.06
HB HGB	g/dl	13.5	13.2			13.5					± 0.2
	g/l	135	132			135					± 2
	mmol/l	8.38	8.22			8.38					± 0.12
HT HCT	%	38.0	38.3			38.0					± 1.0
	l/l	0.380	0.383			0.380					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	238	236			243					± 10
VPM MPV	µm ³ ; fl	8.5	8.2			8.4					± 0.5

ABX Minocal


 (Exp.) **2018-10-05**
 (YYYY - MM - DD)

LOT		CX 426		CAL		(Exp.) 2018-10-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.51								± 0.06
HB HGB	g/dl	13.3	13.3								± 0.2
	g/l	133	133								± 2
	mmol/l	8.26	8.26								± 0.12
HT HCT	%	38.0	38.2								± 1.0
	l/l	0.380	0.382								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	225	235								± 10
VPM MPV	µm ³ ;fl	N/A	8.7								± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A