

# ABX Minocal

**LOT**

**CX 426**

**CAL**



**(Exp.) 2018-10-05**  
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.1</b>	<b>9.1</b>			<b>9.3</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.52</b>	<b>4.51</b>			<b>4.50</b>					± 0.06
HB HGB	g/dl	<b>13.5</b>	<b>13.2</b>			<b>13.5</b>					± 0.2
	g/l	<b>135</b>	<b>132</b>			<b>135</b>					± 2
	mmol/l	<b>8.38</b>	<b>8.22</b>			<b>8.38</b>					± 0.12
HT HCT	%	<b>38.0</b>	<b>38.3</b>			<b>38.0</b>					± 1.0
	l/l	<b>0.380</b>	<b>0.383</b>			<b>0.380</b>					± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>238</b>	<b>236</b>			<b>243</b>					± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>8.5</b>	<b>8.2</b>			<b>8.4</b>					± 0.5

# ABX Minocal

**LOT**

**CX 426**

Rev 1

**CAL**



(Exp.) **2018-10-05**  
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.0</b>	<b>9.0</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.50</b>	<b>4.51</b>							± 0.06
HB HGB	g/dl	<b>13.3</b>	<b>13.3</b>							± 0.2
	g/l	<b>133</b>	<b>133</b>							± 2
	mmol/l	<b>8.26</b>	<b>8.26</b>							± 0.12
HT HCT	%	<b>38.0</b>	<b>38.2</b>							± 1.0
	l/l	<b>0.380</b>	<b>0.382</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>225</b>	<b>235</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>8.7</b>							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A